


Mark[®] 1210c

Where Design Meets Function.

The sleek, elegant design will undoubtedly add style to your practice. Functionality and design engineered into this groundbreaking platform that will add efficiency and profitability to your bottom line.

Engineered and designed with three field replaceable components to *virtually eliminate down time*.

The touch screen graphical interface provides fast, intuitive and flexible workflow. Sound-Eklin's[®] exclusive Sound Vu™ technology results in the best possible images.



The 12" x 10" Cesium detector panel delivers better images at lower technique.

YOUR IMAGE IS OUR BUSINESS.



a VCA ANTECH company

800.268.5354 soundeklin.com

Mark[®] 1210c


Specifications

Sensor Panel

Model	VaTech 1210PCB / Samsung
Type	Amorphous Silicon TFT X-Ray Converter: Csl. TI
Acquisition Time	5 seconds
Imaging Area	10" x 12", actual size 25.4 cm X 30.48 cm, actual size
Pixel Pitch	127 Micron
Pixel Array	2080 x 2560 pixels
Grayscale	14 bit (4,096 levels)
Weight	8.1 lbs 3.7 kg
Dimensions	14.5" w x 15" h x 0.9" d 36.83 cm w x 38.1 cm h x 2.3 cm d

System

Software	eSeries Touch
Weight	41 lbs 18.4 kg
Processor	2.53 GHz Intel Core 2 Duo Processor
Memory	4 GB RAM, 500GB HD
OS	Windows [®] 7
Hard Drive	500GB
DB Backup	16GB
Weight	41 lbs, 18.6 kg
Display	17" Touchscreen



YOUR IMAGE IS OUR BUSINESS.



a VCA ANTECH company

For information or a demonstration call

800.268.5354 soundeklin.com



*The Mark[®]1210c
with Field Replaceable
Electronics Box Removed*