



KONICA MINOLTA

# AeroDR

## VETERINARY SOLUTION

*Speed up and simplify your workflow with AeroDR and ImagePilot.*

*AeroDR is Konica Minolta's lightweight, yet extremely robust flat panel detector, used throughout the world to create high quality digital images. At the same time several thousands of veterinarians worldwide have been using the ImagePilot all-in-one system every day to their full satisfaction. These two tried and trusted products have now been combined to create the AeroDR Veterinary Solution. The ultimate way to speed up and simplify your workflow.*

### Console, Viewer and Archive in one

AeroDR Veterinary Solution offers an easy to use all-in-one Console, Viewer and Archive for your daily veterinary use. With our unique one click acquisition workflow, you simply make the exposure and your image will be displayed almost immediately. Konica Minolta's unique auto processing function AutoPilot guarantees a perfect image every time. With storage capacities up to 2TB, AeroDR Veterinary Solution can archive up to 300,000 images on its expandable server.



### Works with your existing X-ray equipment

Immediately experience the benefits of a faster and more efficient workflow. By using our unique AeroSync technology, your current X-ray system can easily be transformed to a completely wireless DR solution within a matter of minutes. In addition, the AeroDR panel can easily and efficiently be shared between the X-ray room in your practice and a mobile X-ray device. Make the most of your existing equipment and go digital with the AeroDR Veterinary Solution.



Giving Shape to Ideas

# AeroDR

## VETERINARY SOLUTION

### AeroDR Specifications

- CsI Scintillator - High image quality
- Light weight - only 2.9 Kg (35x43 cm), 1.7 Kg (25x30 cm\*)
- Pixel Size - 175µm
- Lithium Ion Capacitor - Fast charging and long life time
- Charging time empty to full - 30 minutes with battery charger
- Carbon Fiber monocoque housing - Robust and durable
- Preview time - 2 seconds



\* available end 2013

### ImagePilot Specifications

- ImagePilot comes equipped with a full set of tools: all you need for your daily practise and more
- Register up to 20 clients and up to 5 Mobile iPad clients
- Standard 380 GB storage ( $\pm$  39.000 images at 35x43 cm), expandable
- AutoPilot automatically applies the correct processing to optimise the image, creating a one click acquisition workflow
- Processing optimised for AeroDR



KONICA MINOLTA

KONICA MINOLTA MEDICAL & GRAPHIC IMAGING EUROPE B.V.  
Frankfurtstraat 40, 1175 RH Lijnden, The Netherlands

Distributed by :