

ZooMax GOLD ELEV

- Veterinary system with four-way elevating table and fixed floor-to-wall tube stand.
- Easy patient positioning, table top height adjustable between 38 and 90 cm. Max. load 100 kg.
- Four-way table top ensures trouble-free operation on the full length of the table.
- Electromagnetic brakes and powerful actuator guarantee safe and convenient work on the table.
- Control of the table with 3 foot switches in one compartment, Up, Down and 4-way top release.

The ZooMax GOLD ELEV system can be supplied with a monoblock generator (3,2-5 kW), or with a conventional high frequency generator (32, 40 or 50 kW).

Please see the price sheet for more information!

