

Veterinary X-Ray Systems for today and tomorrow

ZooMax^{RED}

The **ZooMax RED** veterinary radiographic systems are best suited for veterinary hospitals and clinics. The new approach in administering to pets in a hospital environment makes the task easy, effortless and hygienic for both patient and practitioner. The system allows full access to the patient for easy positioning, restraint and exposure control from either side of the table in trauma situations. This new approach revolutionizes the process of caring for animals.

Fixed SID, but adjustable at install. Longitudinal displacement of the tube is approx. 150 cm. Tube head rotating for oblique studies. Elevating table with 4-way, table top.

