



# ZooMax DIGI GOLD

## Digital Veterinary Radiographic System



All digital with  
**SolidDR™**




Optional VetElev table

- High Power 32kW Digital X-ray System. 220V or 380V connection to mains
- Very short exposure times & dual focus x-ray tube for sharp images
- PC touch console with full APR functions
- Rotatable LED collimator for bright light field
- Electrical 4-way floating table top 152x71cm, cross +/- 10cm, side to side +/- 35cm
- Fastening rails & cleat rope ties for fixing the animal
- Easy withdrawable grid function with storage
- Remote support via internet TeamViewer
- Perfect all digital with Medivet SolidDR
- Quality that will last a long time. Made in Europe

mm	inch
2225	87.6
1885	74.2
1525	60.0
1400	55.1
1000	39.4
820	32.9
735	28.9
730	28.7
560	22.0
550	21.7
420	16.5
350	13.8
245	9.64
120	4.72
100	3.94



Designer	L. Adrovitz	Title <b>ZMX<sup>GOLD</sup> ROOM LAYOUT</b>	Doc. Number		
Drawn by	J. Dudas				
Checked by	L. Breda				
Approved by	L. Adrovitz	File: ZooMax Gold RL-b	Date: 22.05.2012.	Scale: 1:20	Sheet 1 of 1